



User Manual

Digital Duplex PCR dWheSp

Principle

This method describes a routine procedure for the quantitative detection of spelt and wheat DNA. Using the ratio of the concentrations of Positive Droplets (copies/ μ l) together with a conversion factor, the weight of the sample can be determined in percent by weight.



Top population :Spelt
Middle population :Wheat
Bottom population: empty

Contents and Storage

4 tubes of primer mix, lyophilized, for 4x24 reactions. Shipped at ambient temperature, store at -20°C

Supplied on demand: EXCEL table for results calculation

Reagents to be Supplied by User

QX200 ddPCR Supermix for Probes 5ml (Cat No. 186-3024) and other consumables for ddPCR with the QX200 system.

Protocol

1. Add 153 μ l PCR grade water per tube of primer-mix, vortex vigorously and incubate for 5 min at 60°C (store solution at 4°C , stable for 1 week).

2. Add 282 μ l Supermix for probes (2x) and mix well.

Yields 435 μ l ready-to-use Mastermix

3. Mix 17 μ l ready-to-use Mastermix with 5 μ l sample solution in an appropriate PCR reaction vessel. Recommended amount of DNA: 0.6 ng, measured by absorbance and/or 1:32 diluted isolated DNA.

4. Use the following thermal cycling profile:

- 1 10 min, 95°C
- 2 30 s, 95°C
- 3 60 s, 62°C
- 4 Repeat steps 2 to 3 **50 times in total**
- 5 10 min, 98°C
- 6 Hold, 4°C

6. Set the Probes parameter for the droplet reading.

7. Analyse the results with QuantaSoft software. To obtain quantifiable results the droplets/ μ l value for the total amount of DNA (the lowest threshold) should be over 1cp/ μ l. The value of the empty droplets should be under 1cp/ μ l.

Contact

Kantonales Labor Zürich
rene.koepfel@klzh.ch

Further Information

<https://www.microsynth.com/food-testing-assays.html>

Microsynth AG
Schützenstrasse 15
CH-9436 Balgach
Telefon +41 71 722 83 33
Fax +41 71 722 87 58
www.microsynth.com